

MEASURING GEOMETRIES



R-L/ISO 2555

Measuring spindles
in stainless steel 316 L*



L-1 spindle



L-2 spindle



L-3 spindle



L-4 spindle



Axis R-7



R-1 disc



R-2 disc



R-3 disc



R-4 disc



R-5 disc



R-6 disc

* These measuring spindles are intended for measuring viscosity in a 600 mL beaker.
This system allows measurements according to the standard ASTM/ISO 2555.


MEASURING GEOMETRIES



R-L/ISO 2555

Measuring spindles
in stainless steel 316 L*

NAME	REF	DIM. (mm)	VOL (mL)	VISCOSITY RANGE
L-1 spindle	111010	Ø 18.80 - L 65.1	600	15 - 20,000 mPa·s
L-2 spindle	111011	Ø 18.72 - L 6.86	600	50 - 100,000 mPa·s
L-3 spindle	111012	Ø 12.60 - L 1.78	600	200 - 400,000 mPa·s
L-4 spindle	111013	Ø 3.20 - L 31	600	1000 - 2,000,000 mPa·s
Axis R 1-6 without disc	111000	Threaded axis	-	-
R-1 Disc	111001	Ø 56.26	600	3 - 450,000 mPa·s
R-2 Disc	111002	Ø 46.93	600	15 - 1,750,000 mPa·s
R-3 Disc	111003	Ø 34.69	600	30 - 4,400,000 mPa·s
R-4 Disc	111004	Ø 27.30	600	60 - 8,800,000 mPa·s
R-5 Disc	111005	Ø 21.14	600	120 - 17,600,000 mPa·s
R-6 Disc	111006	Ø 14.62	600	260 - 42,600,000 mPa·s
Axis R-7	111007	Ø 3.20	600	1000 - 156,000,000 mPa·s
Axis L-R	111008	Adaptation axis	-	-
"VANE" Measuring spindles				
Vane 72	120017	Ø 21.67 - L 43.38	-	50 - 500,000 mPa·s
Vane 73	111108	Ø 12.67 - L 25.35	-	250 - 25,000,000 mPa·s
Vane 74	111115	Ø 5.89 - L 11.76	-	2,500 - 250,000,000 mPa·s
MK 6 Blades	111105	W 22	-	100 - 20,000,000 mPa·s
"KREBBS" Measuring spindles				
MK-KU 1-10	111100	W 53.98	250	20 - 5,000 mPa·s 40 - 140 KU
MK-75V	111103	W 42.88	250	100 - 50,000 mPa·s
"FANN R1B1" Measuring spindles				
MK-R1B1	119001	Ø 34.49	-	2 - 800,000 mPa·s
MB-R1B1	119002	Ø 36.80	20	-



* These measuring spindles are intended for measuring viscosity in a 600 mL beaker.
This system allows measurements according to the standard ASTM/ISO 2555.