

DIGITAL HEATING MAGNETIC STIRRER



HS-380D Stirrer

Voltage[VAC]	110/220
Frequency[Hz]	50/60
Power[W]	500
Max. stirring quantity(H ₂ O)[L]	20
Stirring point position quantity[mm]	80
Motor type	brushless motor
Speed range[rpm]	200-1500
Speed and temperature display	LCD
Speed display accuracy[rpm]	1
Plate material:	Ceramic coating surface/Stainless steel
Dimension of Plate[mm]	φ155
Heating Temperature Range[°C]	RT-380
Temperature display accuracy[°C]	0.1
Temperature sensor of extra position	PT100(Optional)
Dimension[mm]	250×150×110
Weight[kg]	2.7

Specifications

Model: S-380D

Speed range[rpm]: 200-1500

Speed display accuracy[rpm]: 1

Model: H-380D

Heating Temperature Range[°C]: RT-380

Temperature display accuracy[°C]: 0.1