

## Ultrasonic Cleaner

**Model: ASD-03L**

**Fluxanar Brand / Germany**



### Feature

- All the buttons are touch-sensitive
- Temperature control knob, Max. temperature set up 80°C (cleaning liquid temperature setting), Timer control knob (cleaning time setting), Max. timer set up 30 minutes
- START/PAUSE & Memory button, On/Off button.
- Degas & Sweep button
  - Degassing is the removal of dissolved oxygen trapped in a liquid. In ultrasonic cleaning systems, trapped, dissolved oxygen in the ultrasonic cleaning solution will hinder the ability of the system to clean. This dissolved oxygen must be removed from the system by de-gassing.
  - Sweep function guarantees an optimum sound field distribution in the cleaning bath due to a modulation of the frequencies.

### Specifications

Model	ASD-03L
Ultrasonic power	160 W
Heating power	150W
Ultrasonic frequency	37KHz
Drain	NO
Internal size(mm) L×W×D	238×138×100
Tank volume	3L
Transducer PC	2 x 80w
Unit Weight	3.2 KG