

HIQLAB



Model	FHD-1200D	FHD-1500D	FHD-1800D	FHD-2000D
FA (Fiberglass) Series :	Chemical Resistant Fiber Glass Fume Hood With digital LCD Display			
Model	FHD-1200A	FHD-1500A	FHD-1800A	FHD-2000A
FA Analog Control Series:	Chemical Resistant Fiber Glass Fume Hood With Analog Control			
Work Top:	FRP (Chemical resistant Fiber glass reinforced plastics Work Top/ Phenolic Work Top(Optional))			
Material Inner :	FRP (Chemical resistant Fiber glass reinforced plastics materials)			
Material Outer :	FRP (Chemical resistant Fiber glass reinforced plastics materials)			
Dimensions : (W x D x H mm)	1200 x 900 x 2350	1500 x 900 x 2350	1800 x 900 x 2350	2000 x 900 x 2350
	(โครงสร้างทากจากไฟเบอร์กลาสเสริมแรง FRP ชนิดที่ทนทานต่อการกัดกร่อนสูง(ISO TYPE)			
Work Top :	FRP (Chemical resistant Fiber glass reinforced plastics materials)			
Blower (Sirrcco Fan)	1 HP	1 HP	1 HP	2 HP
Utility Device	Water Cock, Cup Sink Drain, Fluorescent Lamp			
Slide Door	Slide Type Door With Tempered Safety Glass			

