

A200




### Introduction

A200 digital thermo-hygrometer is a kind of measuring instrument, which is developed by HUATO company according to advanced techniques from other countries and customers' requirements. It has a large LCD display, which can display temperature, humidity, date, clock as well as environment status simultaneously. Besides, it has a low battery consumption (1 block 9V battery can last 12 months working time) so that you do not need to exchange battery frequently.

### Features:

- Both °C and °F are available.
- Display temperature, humidity, data, time as well as environment status simultaneously.
- Reliable and stable with high quality and good appearance.
- MAX/MIN function key to show the maximum minimum value.
- Whole point buzzer time function and alarm clock.
- Very easy to use and competitive.
- Desk-top and wall-mounted.

### Standard accessories



### Technical Specifications

Temperature resolution	0.1 °C (0.2°F)
Temperature accuracy	1 °C
Humidity accuracy	5 %RH
Temperature measuring range	-20~70°C / -4~158°F
Humidity measuring range	10~99%RH
Power supply	9V Block battery
Product dimension	135×124×35mm (5.31"×4.48"×1.37" inch)
LCD size	97×78mm (3.81"×1.37" inch)
Weight	330g
Package dimension	152×132×45mm(5.99"×5.2"×1.77" inch)
Accessories	9V block battery, screw, user manual and carrying box .

### Applications

It has been widely applied in medical field, warehouse, office, lab and agricultural industry etc.