

Hot Plate & Stirrer

V-TECH

Model: HS-370



Features

- Chemical resistant ceramic coated
- Constructed for a long & durable life
- Charming feature which statute indicator
- Uniform heating conditions
- Very quick obtainment of desired Temperature.
- PT100 Temperature sensor(Optional)
- Stable stand(Optional)

MODEL	HS-370
Speed Range (RPM)	100-1600
Temperature Range	Ambient +370°C
Temperature controller	Proportional controller
Heater	Interchangeable Heater
Plate area	190 x 190mm
Dimension (W x D x H) (mm)	335 x 230 x 90mm
Power	220V, 50/60Hz