

## Auto-Dry Cabinet

### Dry Active

### Korea Brand



#### Features:

- Stores various products that are sensitive to humidity changes: Cameras, Electronic goods, Reagents, Specimens, Medicines, Sem-conductors etc.
- High-density foam gaskets perfectly prevent any contact with the outside.
- It is lighter than those of other companies, and excellent for protecting fragileness as all windows are made of PC.
- A wide digital LCD display screen indicates internal humidity and temperature.
- Auto-Dry Cabinet devices reliably keep internal humidity at 20%-40%.

#### Technical data:

|                          |   |
|--------------------------|---|
| Model                    | KA.33-73  |
| Size                     | 47.6L   |
| Outside Dimensions       | 300(W)×380(D)×635(H)mm                                    |
| Internal Dimensions      | 280(W)×340(D)×500(H)mm                                    |
| Size of Shelf            | 280(W)×340(D)mm   |
| Accessories              | 3 Shelves, 1 Digital Thermo-Hygrometer(built-in) included |
| Material                 | ABS+PC  |
| Humidity Control Range   | 20-40%  |
| Indicating humidity      | 10-80%  |
| Deviation of humidity    | ±3%   |
| Indicating temperature   | -50°C~99°C(-58°F~210°F)                                   |
| Deviation of temperature | ±1°C(±2°F)  |