

# LABORHITS

## Microliter Pipettes LLG-proMLP

**Precise and Reliable Pipette**

**Extremely light weight while the robustness still remains**


- High accuracy, reliability and reproducibility
- Perfect ergonomics for highest comfort
- Colour-coded push-button for quick operation
- Extremely low pipetting forces
- Built-in volume locking for error-free results
- Tip ejector button\*, adjustable for right- and left-handed users
- Completely autoclavable at 121 °C
- Including individual calibration certificate (according to EN ISO 8655)



\*not available for 5 ml and 10 ml pipettes

Capacity µl	Increments µl	Accuracy max. vol. (≤ ± R%)	Precision max. vol. (≤ VK%)	Colour code	PK	Cat No.
0,1 ... 2	0,002	1,65	0,75	dark orange	1	<b>6.270 420</b>
0,5 ... 10	0,02	1,10	0,44	red	1	<b>6.270 421</b>
2 ... 20	0,02	1,00	0,33	light yellow	1	<b>6.270 422</b>
10 ... 100	0,2	0,80	0,17	light orange	1	<b>6.270 423</b>
20 ... 200	0,2	0,80	0,17	dark yellow	1	<b>6.270 424</b>
100 ... 1000	2	0,80	0,17	dark blue	1	<b>6.270 425</b>
500 ... 5000	2	0,66	0,18	purple	1	<b>6.270 426</b>
100 ... 10000	20	0,60	0,18	light blue	1	<b>6.270 427</b>

**NEW**



### LLG-Filter tips Low Retention 2.0, PP

Available in easy-to-open boxes with color coded racks. Compatible with all common pipette brands. Improved optical clarity and pipette fit. Colored coded rack inserts for easy identification. Graduated, sterile, DNase and RNase free

Capacity µl	Package contents	PK	Cat. No.
0,1 ... 10	10 x 96	960	<b>4.668 784</b>
0,1 ... 10*	10 x 96	960	<b>4.668 785</b>
1 ... 20	10 x 96	960	<b>4.668 786</b>
1 ... 100	10 x 96	960	<b>4.668 787</b>
1 ... 200	10 x 96	960	<b>4.668 788</b>
100 ... 1000	8 x 96	768	<b>4.668 789</b>

\* Tips XL