

Viscometer



Model	VT-0100K	VT-2000K
Measurement Tolerance	± 2% (Newtonian liquid)	
Rotor Specification	1, 2, 3, 4# (0# optional)	
Measurement Range	1-100,000 mpa.s	1-2,000,000 mpa.s
Rotating Speed of Rotor	6 / 12 / 30 / 60 rpm	0.3 / 0.6 / 1.5 / 3 / 6 / 12 / 30 / 60 rpm
Material breakdown	ABS HOUSE, STAINLESS IRON, STEEL	
Power	AC (220-240VAC 50Hz)	
Working Condition	Temperature 5-35°C Relative humidity no more than 80%	
Interface	No	
Outline dimension (L×W×H)	14.6"×12.8"×11" (370×325×280 mm)	
Automatic Switch	free selection of proper rotating speed or rotor number	
Stable Cursor	Strip-type cursor fully charged displayed when reading basically stable	

Distributor by Siripanya Group

Tel.086-6540217, 088-2585904, 091-8515297

E-mail: Admin1@panyalab.com, Admin2@panyalab.com,

puripim@hotmail.com

www.panyalab.com