

- Voltage: 110-240 (VAC)
- Speed: 2200RPM
- Power: 6.6W



Magnetic Stirrer Model: SI-2200

- Wide voltage input 110-240 (VAC), high voltage adaptability, energy saving and durable.
- Rotating speed: 0-2200RPM.
- Fully enclosed flame retardant plastic enclosure.
- Plastic tray can be anti-splash water and anti-rust.

- Rotating Speed: 2850RPM
- Power: 60W



Vortex Mixer Model: VT-60

- Suitable for short time operation or long time continuous work.
- Running speed, 2850RPM.
- Compact, super shockproof, suitable for high speed operation.
- Stable and reliable