



- Voltage: 100-240 (VAC)
- Frequency: 50/60Hz
- Power: 20W
- Turnover diameter: 22mm
- Speed: 10-200RPM
- Maximum load: 10kg
- Operation mode: Orbital

Orbital Shaker Model : OS-200

- Wide voltage input, high adaptability of voltage fluctuation, energy saving and durable.
- LED display of time, speed, temperature.
- Speed adjustable, time controllable.
- Simple design, easy operation, quick reaction.
- Small size, work in humid environment, in the refrigerator and the incubator.
- Design with brush-less motor, long service life.
- Smooth and stable operation.



- Voltage: 100-240 (VAC)
- Frequency: 50/60Hz
- Power: 20W
- Swing: $\pm 9^\circ$
- Speed: 10-100RPM
- Maximum load: 10kg
- Operation mode: Rocker

Rocker shaker Model: RS-100

- Wide voltage input, high adaptability of voltage fluctuation, energy saving and durable.
- LED display of time, speed, temperature.
- Speed adjustable, time controllable.
- Simple design, easy operation, quick reaction.
- Small size, work in humid environment, in the refrigerator and the incubator.
- Design with brush-less motor, long service life.
- Smooth and stable operation.