



PROMOTION

P-CLEAN

- ✓ STERILITY
- ✓ NUCLEASE-FREE
- ✓ PYROGEN-FREE



PRODUCT

- Closures without liner, leakproof guarantee
- Undergo pressure of 80 kPa, no leakage
- Produced in cleanroom class ISO 8, Endotoxin less than 0.125 EU/ml
- Chemical resistant
- Torque design on lip, which for auto-filling line

HDPE Wide Mouth Reagent Bottles, PP screw closure

Capability	Height (mm)	O.D. (mm)	Natural 	Packing	ราคาต่อยกแพ็ค	ราคาต่อชิ้น
250 ml	128	60	•	50 pcs/pack	1,850.-	50.-
500 ml	163	72	•	50 pcs/pack	2,750.-	75.-
1000 ml	193	91	•	50 pcs/pack	3,700.-	95.-

Wide Mouth Reagent Bottles, PP screw closure

Capability	Height (mm)	O.D. (mm)	Natural 	Packing	ราคาต่อยกแพ็ค	ราคาต่อชิ้น
125 ml	96	50	•	50 pcs/pack	1,250.-	45.-
250 ml	128	60	•	50 pcs/pack	1,450.-	55.-
500 ml	163	72	•	50 pcs/pack	2,650.-	70.-
1000 ml	193	91	•	50 pcs/pack	3,350.-	85.-

HDPE Wide Mouth Reagent Bottles, PP screw closure

Capability	Height mm	O.D. mm	Natural 	Packing	ราคาต่อยกแพ็ค	ราคาต่อชิ้น
250 ml	128	60	•	50 pcs/pack	2,500.-	58.-
500 ml	169	72	•	50 pcs/pack	3,650.-	83.-
1000 ml	198	90	•	50 pcs/pack	4,350.-	90.-

**สินค้ามีจำนวนจำกัด ราคาดังกล่าวยังไม่รวมภาษีมูลค่าเพิ่ม 7%
 ห้างฯ ขอสงวนสิทธิ์กรณีสั่งซื้อสินค้าไม่ถึง 1,000 บาท คิดค่าขนส่งตามจริง
 ห้างฯ ขอสงวนสิทธิ์ในการเปลี่ยนแปลงโดยไม่ต้องแจ้งให้ทราบล่วงหน้า