

## VACUUM PUMP MODEL TVP-20



- Speed of evacuation ( L / Min ) : 20
- Ultimate pressure vacuum :  $\geq 0.08$  Mpa  
200 mbar
- Voltagerating : 220V / 110V; 50Hz / 60Hz
- Motor power (w) : 160W
- Inlet (mm) :  $\varnothing 6$
- Outlet (mm) : InnerSilencer
- Temperature of working environment ( °C ) : 7 - 40
- Function : Vacuum
- Pumphead : 1
- Temp of the body ( °C ) : < 55
- Dimensions ( LXBXH ) (mm) : 215 x 120 x 235
- Noiselevel ( D B ) : < 60
- Weight ( kg ) : 7.5
- Diaphragm : NBR
- Valves : NBR