



IMPACT IMPACT IMPACT

ISO 3310

Woven Wire Mesh

ASTM E11

Half Height Frames

BS 410

Perforated Plate

ISO 5223

Guaranteed to meet International Standards

75mm

100mm

150mm

200mm

250mm

300mm

315mm

350mm

400mm

450mm

8 inch

12 inch

Test Sieves

Index

| Item | Page |
|--|--------------|
| American Standard Sieves to ASTM E11-04 | 7 |
| American Standard Sieves (Half-Height) | 8 |
| European Standard Sieves to ISO 3310/1 and BS 410/1 | 3 - 4 |
| European Standard Sieves to ISO 3310/2 and BS 410/2 | 4 - 5 |
| European Standard Sieves (Half-Height) | 6 - 7 |
| Grain Sieves to ISO 5223 | 8 |
| Large Capacity Screen Tray Sieves | 9 |
| Sieve Specifications | 10 |
| Wet Washing Sieves to ISO 3310/1 and ASTM E11-04 | 10 |

European Standard Sieves to ISO 3310/1 and BS 410/1

Woven Wire Mesh

| Aperture | Frame Diameter | | | | | | | | | |
|-------------------|----------------|-------|-------|-------|-------|-------|-------|-------|-------|-------|
| | 75mm | 100mm | 150mm | 200mm | 250mm | 300mm | 315mm | 350mm | 400mm | 450mm |
| 20µm woven wire | SX941 | SX100 | SX001 | SV496 | SX200 | SV596 | SX300 | SX400 | SX500 | SV280 |
| 25µm woven wire | SX942 | SX101 | SX002 | SV497 | SX201 | SV597 | SX301 | SX401 | SX501 | SV281 |
| 32µm woven wire | SX943 | SX102 | SX003 | SV498 | SX202 | SV598 | SX302 | SX402 | SX502 | SV282 |
| 36µm woven wire | SX944 | SX103 | SX004 | SV499 | SX203 | SV599 | SX303 | SX403 | SX503 | SV343 |
| 38µm woven wire | SX945 | SX104 | SX005 | SV500 | SX204 | SV600 | SX304 | SX404 | SX504 | SV283 |
| 40µm woven wire | SX946 | SX105 | SX006 | SV501 | SX205 | SV601 | SX305 | SX405 | SX505 | SV344 |
| 45µm woven wire | SX947 | SX106 | SX007 | SV502 | SX206 | SV602 | SX306 | SX406 | SX506 | SV284 |
| 50µm woven wire | SX948 | SX107 | SX008 | SV503 | SX207 | SV603 | SX307 | SX407 | SX507 | SV345 |
| 53µm woven wire | SX949 | SX108 | SX009 | SV504 | SX208 | SV604 | SX308 | SX408 | SX508 | SV285 |
| 56µm woven wire | SX950 | SX109 | SX010 | SV505 | SX209 | SV605 | SX309 | SX409 | SX509 | SV346 |
| 63µm woven wire | SX951 | SX110 | SX011 | SV506 | SX210 | SV606 | SX310 | SX410 | SX510 | SV286 |
| 71µm woven wire | SX952 | SX111 | SX012 | SV507 | SX211 | SV607 | SX311 | SX411 | SX511 | SV347 |
| 75µm woven wire | SX953 | SX112 | SX013 | SV508 | SX212 | SV608 | SX312 | SX412 | SX512 | SV287 |
| 80µm woven wire | SX954 | SX113 | SX014 | SV509 | SX213 | SV609 | SX313 | SX413 | SX513 | SV348 |
| 90µm woven wire | SX955 | SX114 | SX015 | SV510 | SX214 | SV610 | SX314 | SX414 | SX514 | SV288 |
| 100µm woven wire | SX956 | SX115 | SX016 | SV511 | SX215 | SV611 | SX315 | SX415 | SX515 | SV349 |
| 106µm woven wire | SX957 | SX116 | SX017 | SV512 | SX216 | SV612 | SX316 | SX416 | SX516 | SV289 |
| 112µm woven wire | SX958 | SX117 | SX018 | SV513 | SX217 | SV613 | SX317 | SX417 | SX517 | SV350 |
| 125µm woven wire | SX959 | SX118 | SX019 | SV514 | SX218 | SV614 | SX318 | SX418 | SX518 | SV290 |
| 140µm woven wire | SX960 | SX119 | SX020 | SV515 | SX219 | SV615 | SX319 | SX419 | SX519 | SV351 |
| 150µm woven wire | SX961 | SX120 | SX021 | SV516 | SX220 | SV616 | SX320 | SX420 | SX520 | SV291 |
| 160µm woven wire | SX962 | SX121 | SX022 | SV517 | SX221 | SV617 | SX321 | SX421 | SX521 | SV352 |
| 180µm woven wire | SX963 | SX122 | SX023 | SV518 | SX222 | SV618 | SX322 | SX422 | SX522 | SV292 |
| 200µm woven wire | SX964 | SX123 | SX024 | SV519 | SX223 | SV619 | SX323 | SX423 | SX523 | SV353 |
| 212µm woven wire | SX965 | SX124 | SX025 | SV520 | SX224 | SV620 | SX324 | SX424 | SX524 | SV293 |
| 224µm woven wire | SX966 | SX125 | SX026 | SV521 | SX225 | SV621 | SX325 | SX425 | SX525 | SV354 |
| 250µm woven wire | SX967 | SX126 | SX027 | SV522 | SX226 | SV622 | SX326 | SX426 | SX526 | SV294 |
| 280µm woven wire | SX968 | SX127 | SX028 | SV523 | SX227 | SV623 | SX327 | SX427 | SX527 | SV355 |
| 300µm woven wire | SX969 | SX128 | SX029 | SV524 | SX228 | SV624 | SX328 | SX428 | SX528 | SV295 |
| 315µm woven wire | SX970 | SX129 | SX030 | SV525 | SX229 | SV625 | SX329 | SX429 | SX529 | SV356 |
| 355µm woven wire | SX971 | SX130 | SX031 | SV526 | SX230 | SV626 | SX330 | SX430 | SX530 | SV296 |
| 400µm woven wire | SX972 | SX131 | SX032 | SV527 | SX231 | SV627 | SX331 | SX431 | SX531 | SV357 |
| 425µm woven wire | SX973 | SX132 | SX033 | SV528 | SX232 | SV628 | SX332 | SX432 | SX532 | SV297 |
| 450µm woven wire | SX974 | SX133 | SX034 | SV529 | SX233 | SV629 | SX333 | SX433 | SX533 | SV358 |
| 500µm woven wire | SX975 | SX134 | SX035 | SV530 | SX234 | SV630 | SX334 | SX434 | SX534 | SV298 |
| 560µm woven wire | SX976 | SX135 | SX036 | SV531 | SX235 | SV631 | SX335 | SX435 | SX535 | SV359 |
| 600µm woven wire | SX977 | SX136 | SX037 | SV532 | SX236 | SV632 | SX336 | SX436 | SX536 | SV299 |
| 630µm woven wire | SX978 | SX137 | SX038 | SV533 | SX237 | SV633 | SX337 | SX437 | SX537 | SV360 |
| 710µm woven wire | SX979 | SX138 | SX039 | SV534 | SX238 | SV634 | SX338 | SX438 | SX538 | SV300 |
| 800µm woven wire | SX980 | SX139 | SX040 | SV535 | SX239 | SV635 | SX339 | SX439 | SX539 | SV361 |
| 850µm woven wire | SX981 | SX140 | SX041 | SV536 | SX240 | SV636 | SX340 | SX440 | SX540 | SV301 |
| 900µm woven wire | SX982 | SX141 | SX042 | SV537 | SX241 | SV637 | SX341 | SX441 | SX541 | SV362 |
| 1.00mm woven wire | SX983 | SX142 | SX043 | SV538 | SX242 | SV638 | SX342 | SX442 | SX542 | SV302 |
| 1.12mm woven wire | SX984 | SX143 | SX044 | SV539 | SX243 | SV639 | SX343 | SX443 | SX543 | SV363 |
| 1.18mm woven wire | SX985 | SX144 | SX045 | SV540 | SX244 | SV640 | SX344 | SX444 | SX544 | SV303 |
| 1.25mm woven wire | SX986 | SX145 | SX046 | SV541 | SX245 | SV641 | SX345 | SX445 | SX545 | SV364 |
| 1.40mm woven wire | SX987 | SX146 | SX047 | SV542 | SX246 | SV642 | SX346 | SX446 | SX546 | SV304 |
| 1.60mm woven wire | SX988 | SX147 | SX048 | SV543 | SX247 | SV643 | SX347 | SX447 | SX547 | SV365 |
| 1.70mm woven wire | SX989 | SX148 | SX049 | SV544 | SX248 | SV644 | SX348 | SX448 | SX548 | SV305 |
| 1.80mm woven wire | SX990 | SX149 | SX050 | SV545 | SX249 | SV645 | SX349 | SX449 | SX549 | SV366 |
| 2.00mm woven wire | SX991 | SX150 | SX051 | SV546 | SX250 | SV646 | SX350 | SX450 | SX550 | SV306 |
| 2.24mm woven wire | SX992 | SX151 | SX052 | SV547 | SX251 | SV647 | SX351 | SX451 | SX551 | SV367 |
| 2.36mm woven wire | SX993 | SX152 | SX053 | SV548 | SX252 | SV648 | SX352 | SX452 | SX552 | SV307 |
| 2.50mm woven wire | SX994 | SX153 | SX054 | SV549 | SX253 | SV649 | SX353 | SX453 | SX553 | SV368 |
| 2.80mm woven wire | SX995 | SX154 | SX055 | SV550 | SX254 | SV650 | SX354 | SX454 | SX554 | SV308 |
| 3.15mm woven wire | | | SX056 | SV551 | SX255 | SV651 | SX355 | SX455 | SX555 | SV369 |
| 3.35mm woven wire | | | SX057 | SV552 | SX256 | SV652 | SX356 | SX456 | SX556 | SV309 |
| 3.55mm woven wire | | | SX058 | SV553 | SX257 | SV653 | SX357 | SX457 | SX557 | SV370 |
| 4.00mm woven wire | | | SX059 | SV554 | SX258 | SV654 | SX358 | SX458 | SX558 | SV310 |
| 4.50mm woven wire | | | SX060 | SV555 | SX259 | SV655 | SX359 | SX459 | SX559 | SV371 |

European Standard Sieves to ISO 3310/1 and BS 410/1

Woven Wire Mesh

| Aperture | Frame Diameter | | | | | | | |
|--------------------|----------------|-------|-------|-------|-------|-------|-------|-------|
| | 150mm | 200mm | 250mm | 300mm | 315mm | 350mm | 400mm | 450mm |
| 4.75mm woven wire | SX061 | SV556 | SX260 | SV656 | SX360 | SX460 | SX560 | SV311 |
| 5.00mm woven wire | SX062 | SV557 | SX261 | SV657 | SX361 | SX461 | SX561 | SV372 |
| 5.60mm woven wire | SX063 | SV558 | SX262 | SV658 | SX362 | SX462 | SX562 | SV312 |
| 6.30mm woven wire | SX064 | SV559 | SX263 | SV659 | SX363 | SX463 | SX563 | SV373 |
| 6.70mm woven wire | SX065 | SV560 | SX264 | SV660 | SX364 | SX464 | SX564 | SV313 |
| 7.10mm woven wire | SX066 | SV561 | SX265 | SV661 | SX365 | SX465 | SX565 | SV374 |
| 8.00mm woven wire | SX067 | SV562 | SX266 | SV662 | SX366 | SX466 | SX566 | SV314 |
| 9.00mm woven wire | SX068 | SV563 | SX267 | SV663 | SX367 | SX467 | SX567 | SV375 |
| 9.50mm woven wire | SX069 | SV564 | SX268 | SV664 | SX368 | SX468 | SX568 | SV315 |
| 10.0mm woven wire | SX070 | SV565 | SX269 | SV665 | SX369 | SX469 | SX569 | SV376 |
| 11.2mm woven wire | SX071 | SV566 | SX270 | SV666 | SX370 | SX470 | SX570 | SV316 |
| 12.5mm woven wire | SX072 | SV567 | SX271 | SV667 | SX371 | SX471 | SX571 | SV377 |
| 13.2mm woven wire | SX073 | SV568 | SX272 | SV668 | SX372 | SX472 | SX572 | SV317 |
| 14.0mm woven wire | SX074 | SV569 | SX273 | SV669 | SX373 | SX473 | SX573 | SV378 |
| 16.0mm woven wire | SX075 | SV570 | SX274 | SV670 | SX374 | SX474 | SX574 | SV318 |
| 18.0mm woven wire | SX076 | SV571 | SX275 | SV671 | SX375 | SX475 | SX575 | SV319 |
| 19.0mm woven wire | SX077 | SV572 | SX276 | SV672 | SX376 | SX476 | SX576 | SV320 |
| 20.0mm woven wire | SX078 | SV573 | SX277 | SV673 | SX377 | SX477 | SX577 | SV321 |
| 22.4mm woven wire | SX079 | SV574 | SX278 | SV674 | SX378 | SX478 | SX578 | SV322 |
| 25.0mm woven wire | SX080 | SV575 | SX279 | SV675 | SX379 | SX479 | SX579 | SV323 |
| 26.5mm woven wire | SX081 | SV576 | SX280 | SV676 | SX380 | SX480 | SX580 | SV324 |
| 28.0mm woven wire | SX082 | SV577 | SX281 | SV677 | SX381 | SX481 | SX581 | SV325 |
| 31.5mm woven wire | SX083 | SV578 | SX282 | SV678 | SX382 | SX482 | SX582 | SV326 |
| 35.5mm woven wire | SX084 | SV579 | SX283 | SV679 | SX383 | SX483 | SX583 | SV327 |
| 37.5mm woven wire | SX085 | SV580 | SX284 | SV680 | SX384 | SX484 | SX584 | SV328 |
| 40.0mm woven wire | SX086 | SV581 | SX285 | SV681 | SX385 | SX485 | SX585 | SV329 |
| 45.0mm woven wire | SX087 | SV582 | SX286 | SV682 | SX386 | SX486 | SX586 | SV330 |
| 50.0mm woven wire | SX088 | SV583 | SX287 | SV683 | SX387 | SX487 | SX587 | SV331 |
| 53.0mm woven wire | SX089 | SV584 | SX288 | SV684 | SX388 | SX488 | SX588 | SV332 |
| 56.0mm woven wire | SX090 | SV585 | SX289 | SV685 | SX389 | SX489 | SX589 | SV333 |
| 63.0mm woven wire | SX091 | SV586 | SX290 | SV686 | SX390 | SX490 | SX590 | SV334 |
| 71.0mm woven wire | SX092 | SV587 | SX291 | SV687 | SX391 | SX491 | SX591 | SV335 |
| 75.0mm woven wire | SX093 | SV588 | SX292 | SV688 | SX392 | SX492 | SX592 | SV336 |
| 80.0mm woven wire | SX094 | SV589 | SX293 | SV689 | SX393 | SX493 | SX593 | SV337 |
| 90.0mm woven wire | SX095 | SV590 | SX294 | SV690 | SX394 | SX494 | SX594 | SV338 |
| 100.0mm woven wire | SX096 | SV591 | SX295 | SV691 | SX395 | SX495 | SX595 | SV339 |
| 106.0mm woven wire | SX097 | SV592 | SX296 | SV692 | SX396 | SX496 | SX596 | SV340 |
| 112.0mm woven wire | SX098 | SV593 | SX297 | SV693 | SX397 | SX497 | SX597 | SV341 |
| 125.0mm woven wire | SX099 | SV594 | SX298 | SV694 | SX398 | SX498 | SX598 | SV342 |

European Standard Sieves to ISO 3310/2 and BS 410/2

Square Hole Perforated Plate

| Aperture | Frame Diameter | | | | | | | |
|-------------------------|----------------|-------|-------|-------|-------|-------|-------|--|
| | 200mm | 250mm | 300mm | 315mm | 350mm | 400mm | 450mm | |
| 4.00mm perforated plate | SV695 | SX600 | SV736 | SX700 | SX800 | SX900 | SV402 | |
| 4.50mm perforated plate | SV696 | SX601 | SV737 | SX701 | SX801 | SX901 | SV403 | |
| 4.75mm perforated plate | SV697 | SX602 | SV738 | SX702 | SX802 | SX902 | SV404 | |
| 5.00mm perforated plate | SV698 | SX603 | SV739 | SX703 | SX803 | SX903 | SV405 | |
| 5.60mm perforated plate | SV699 | SX604 | SV740 | SX704 | SX804 | SX904 | SV406 | |
| 6.30mm perforated plate | SV700 | SX605 | SV741 | SX705 | SX805 | SX905 | SV407 | |
| 6.70mm perforated plate | SV701 | SX606 | SV742 | SX706 | SX806 | SX906 | SV408 | |
| 7.10mm perforated plate | SV702 | SX607 | SV743 | SX707 | SX807 | SX907 | SV409 | |
| 8.00mm perforated plate | SV703 | SX608 | SV744 | SX708 | SX808 | SX908 | SV410 | |
| 9.00mm perforated plate | SV704 | SX609 | SV745 | SX709 | SX809 | SX909 | SV411 | |
| 9.50mm perforated plate | SV705 | SX610 | SV746 | SX710 | SX810 | SX910 | SV412 | |
| 10.0mm perforated plate | SV706 | SX611 | SV747 | SX711 | SX811 | SX911 | SV413 | |
| 11.2mm perforated plate | SV707 | SX612 | SV748 | SX712 | SX812 | SX912 | SV414 | |
| 12.5mm perforated plate | SV708 | SX613 | SV749 | SX713 | SX813 | SX913 | SV415 | |
| 13.2mm perforated plate | SV709 | SX614 | SV750 | SX714 | SX814 | SX914 | SV416 | |

European Standard Sieves to ISO 3310/2 and BS 410/2

Square Hole Perforated Plate

| Aperture | Frame Diameter | | | | | | |
|-------------------------|----------------|-------|-------|-------|-------|-------|-------|
| | 200mm | 250mm | 300mm | 315mm | 350mm | 400mm | 450mm |
| 14.0mm perforated plate | SV710 | SX615 | SV751 | SX715 | SX815 | SX915 | SV417 |
| 16.0mm perforated plate | SV711 | SX616 | SV752 | SX716 | SX816 | SX916 | SV418 |
| 18.0mm perforated plate | SV712 | SX617 | SV753 | SX717 | SX817 | SX917 | SV419 |
| 19.0mm perforated plate | SV713 | SX618 | SV754 | SX718 | SX818 | SX918 | SV420 |
| 20.0mm perforated plate | SV714 | SX619 | SV755 | SX719 | SX819 | SX919 | SV421 |
| 22.4mm perforated plate | SV715 | SX620 | SV756 | SX720 | SX820 | SX920 | SV422 |
| 25.0mm perforated plate | SV716 | SX621 | SV757 | SX721 | SX821 | SX921 | SV423 |
| 26.5mm perforated plate | SV717 | SX622 | SV758 | SX722 | SX822 | SX922 | SV424 |
| 28.0mm perforated plate | SV718 | SX623 | SV759 | SX723 | SX823 | SX923 | SV425 |
| 31.5mm perforated plate | SV719 | SX624 | SV760 | SX724 | SX824 | SX924 | SV426 |
| 35.5mm perforated plate | SV720 | SX625 | SV761 | SX725 | SX825 | SX925 | SV427 |
| 37.5mm perforated plate | SV721 | SX626 | SV762 | SX726 | SX826 | SX926 | SV428 |
| 40.0mm perforated plate | SV722 | SX627 | SV763 | SX727 | SX827 | SX927 | SV429 |
| 45.0mm perforated plate | SV723 | SX628 | SV764 | SX728 | SX828 | SX928 | SV430 |
| 50.0mm perforated plate | SV724 | SX629 | SV765 | SX729 | SX829 | SX929 | SV431 |
| 53.0mm perforated plate | SV725 | SX630 | SV766 | SX730 | SX830 | SX930 | SV432 |
| 56.0mm perforated plate | SV726 | SX631 | SV767 | SX731 | SX831 | SX931 | SV433 |
| 63.0mm perforated plate | SV727 | SX632 | SV768 | SX732 | SX832 | SX932 | SV434 |
| 71.0mm perforated plate | SV728 | SX633 | SV769 | SX733 | SX833 | SX933 | SV435 |
| 75.0mm perforated plate | SV729 | SX634 | SV770 | SX734 | SX834 | SX934 | SV436 |
| 80.0mm perforated plate | SV730 | SX635 | SV771 | SX735 | SX835 | SX935 | SV437 |
| 90.0mm perforated plate | SV731 | SX636 | SV772 | SX736 | SX836 | SX936 | SV438 |
| 100mm perforated plate | SV732 | SX637 | SV773 | SX737 | SX837 | SX937 | SV439 |
| 106mm perforated plate | SV733 | SX638 | SV774 | SX738 | SX838 | SX938 | SV440 |
| 112mm perforated plate | SV734 | SX639 | SV775 | SX739 | SX839 | SX939 | SV441 |
| 125mm perforated plate | SV735 | SX640 | SV776 | SX740 | SX840 | SX940 | SV442 |

Round Hole Perforated Plate

| Aperture | Frame Diameter | |
|------------------------------------|----------------|--------|
| | 200mm | 300mm |
| 4.00mm round hole perforated plate | SV695R | SV736R |
| 4.50mm round hole perforated plate | SV696R | SV737R |
| 4.75mm round hole perforated plate | SV697R | SV738R |
| 5.00mm round hole perforated plate | SV698R | SV739R |
| 5.60mm round hole perforated plate | SV699R | SV740R |
| 6.30mm round hole perforated plate | SV700R | SV741R |
| 6.70mm round hole perforated plate | SV701R | SV742R |
| 7.10mm round hole perforated plate | SV702R | SV743R |
| 8.00mm round hole perforated plate | SV703R | SV744R |
| 9.00mm round hole perforated plate | SV704R | SV745R |
| 9.50mm round hole perforated plate | SV705R | SV746R |
| 10.0mm round hole perforated plate | SV706R | SV747R |
| 11.2mm round hole perforated plate | SV707R | SV748R |
| 12.5mm round hole perforated plate | SV708R | SV749R |
| 13.2mm round hole perforated plate | SV709R | SV750R |
| 14.0mm round hole perforated plate | SV710R | SV751R |
| 16.0mm round hole perforated plate | SV711R | SV752R |
| 18.0mm round hole perforated plate | SV712R | SV753R |
| 19.0mm round hole perforated plate | SV713R | SV754R |
| 20.0mm round hole perforated plate | SV714R | SV755R |
| 22.4mm round hole perforated plate | SV715R | SV756R |

| Aperture | Frame Diameter | |
|------------------------------------|----------------|--------|
| | 200mm | 300mm |
| 25.0mm round hole perforated plate | SV716R | SV757R |
| 26.5mm round hole perforated plate | SV717R | SV758R |
| 28.0mm round hole perforated plate | SV718R | SV759R |
| 31.5mm round hole perforated plate | SV719R | SV760R |
| 35.5mm round hole perforated plate | SV720R | SV761R |
| 37.5mm round hole perforated plate | SV721R | SV762R |
| 40.0mm round hole perforated plate | SV722R | SV763R |
| 45.0mm round hole perforated plate | SV723R | SV764R |
| 50.0mm round hole perforated plate | SV724R | SV765R |
| 53.0mm round hole perforated plate | SV725R | SV766R |
| 56.0mm round hole perforated plate | SV726R | SV767R |
| 63.0mm round hole perforated plate | SV727R | SV768R |
| 71.0mm round hole perforated plate | SV728R | SV769R |
| 75.0mm round hole perforated plate | SV729R | SV770R |
| 80.0mm round hole perforated plate | SV730R | SV771R |
| 90.0mm round hole perforated plate | SV731R | SV772R |
| 100mm round hole perforated plate | SV732R | SV773R |
| 106mm round hole perforated plate | SV733R | SV774R |
| 112mm round hole perforated plate | SV734R | SV775R |
| 125mm round hole perforated plate | SV735R | SV776R |

Lids and Receivers

| | Frame Diameter | | | | | | | | | |
|----------------------------|----------------|-------|-------|-------|-------|-------|-------|-------|-------|-------|
| | 75mm | 100mm | 150mm | 200mm | 250mm | 300mm | 315mm | 350mm | 400mm | 450mm |
| Lid | SV391 | SV467 | SV451 | SV470 | SV479 | SV473 | SV482 | SV485 | SV488 | SV476 |
| Receiver | SV392 | SV468 | SV452 | SV471 | SV480 | SV474 | SV483 | SV486 | SV489 | SV477 |
| Receiver/Separator | SV393 | SV469 | SV453 | SV472 | SV481 | SV475 | SV484 | SV487 | SV490 | SV478 |
| Wet Washing Lid & Receiver | SV390 | SV454 | SV456 | SV461 | SV462 | SV463 | SV464 | SV457 | SV465 | SV466 |



European Standard Sieves to ISO 3310/1 and BS 410/1

Half Height Frames

Woven Wire Mesh

| Aperture | – Frame Diameter – | |
|--------------------|--------------------|--------|
| | 200mm | 300mm |
| 20µm woven wire | HSV496 | HSV596 |
| 25µm woven wire | HSV497 | HSV597 |
| 32µm woven wire | HSV498 | HSV598 |
| 36µm woven wire | HSV499 | HSV599 |
| 38µm woven wire | HSV500 | HSV600 |
| 40µm woven wire | HSV501 | HSV601 |
| 45µm woven wire | HSV502 | HSV602 |
| 50µm woven wire | HSV503 | HSV603 |
| 53µm woven wire | HSV504 | HSV604 |
| 56µm woven wire | HSV505 | HSV605 |
| 63µm woven wire | HSV506 | HSV606 |
| 71µm woven wire | HSV507 | HSV607 |
| 75µm woven wire | HSV508 | HSV608 |
| 80µm woven wire | HSV509 | HSV609 |
| 90µm woven wire | HSV510 | HSV610 |
| 100µm woven wire | HSV511 | HSV611 |
| 106µm woven wire | HSV512 | HSV612 |
| 112µm woven wire | HSV513 | HSV613 |
| 125µm woven wire | HSV514 | HSV614 |
| 140µm woven wire | HSV515 | HSV615 |
| 150µm woven wire | HSV516 | HSV616 |
| 160µm woven wire | HSV517 | HSV617 |
| 180µm woven wire | HSV518 | HSV618 |
| 200µm woven wire | HSV519 | HSV619 |
| 212µm woven wire | HSV520 | HSV620 |
| 224µm woven wire | HSV521 | HSV621 |
| 250µm woven wire | HSV522 | HSV622 |
| 280µm woven wire | HSV523 | HSV623 |
| 300µm woven wire | HSV524 | HSV624 |
| 315µm woven wire | HSV525 | HSV625 |
| 355µm woven wire | HSV526 | HSV626 |
| 400µm woven wire | HSV527 | HSV627 |
| 425µm woven wire | HSV528 | HSV628 |
| 450µm woven wire | HSV529 | HSV629 |
| 500µm woven wire | HSV530 | HSV630 |
| 560µm woven wire | HSV531 | HSV631 |
| 600µm woven wire | HSV532 | HSV632 |
| 630µm woven wire | HSV533 | HSV633 |
| 710µm woven wire | HSV534 | HSV634 |
| 800µm woven wire | HSV535 | HSV635 |
| 850µm woven wire | HSV536 | HSV636 |
| 900µm woven wire | HSV537 | HSV637 |
| 1.00mm woven wire | HSV538 | HSV638 |
| 1.12mm woven wire | HSV539 | HSV639 |
| 1.18mm woven wire | HSV540 | HSV640 |
| 1.25mm woven wire | HSV541 | HSV641 |
| 1.4mm woven wire | HSV542 | HSV642 |
| 1.60 mm woven wire | HSV543 | HSV643 |
| 1.70 mm woven wire | HSV544 | HSV644 |
| 1.80mm woven wire | HSV545 | HSV645 |

| Aperture | – Frame Diameter – | |
|---------------------|--------------------|--------|
| | 200mm | 300mm |
| 2.00 mm woven wire | HSV546 | HSV646 |
| 2.24 mm woven wire | HSV547 | HSV647 |
| 2.36mm woven wire | HSV548 | HSV648 |
| 2.50 mm woven wire | HSV549 | HSV649 |
| 2.80 mm woven wire | HSV550 | HSV650 |
| 3.15 mm woven wire | HSV551 | HSV651 |
| 3.35 mm woven wire | HSV552 | HSV652 |
| 3.55mm woven wire | HSV553 | HSV653 |
| 4.00 mm woven wire | HSV554 | HSV654 |
| 4.50 mm woven wire | HSV555 | HSV655 |
| 4.75mm woven wire | HSV556 | HSV656 |
| 5.00 mm woven wire | HSV557 | HSV657 |
| 5.60 mm woven wire | HSV558 | HSV658 |
| 6.30 mm woven wire | HSV559 | HSV659 |
| 6.70 mm woven wire | HSV560 | HSV660 |
| 7.10 mm woven wire | HSV561 | HSV661 |
| 8.00mm woven wire | HSV562 | HSV662 |
| 9.00 mm woven wire | HSV563 | HSV663 |
| 9.50 mm woven wire | HSV564 | HSV664 |
| 10.0 mm woven wire | HSV565 | HSV665 |
| 11.2mm woven wire | HSV566 | HSV666 |
| 12.5mm woven wire | HSV567 | HSV667 |
| 13.2 mm woven wire | HSV568 | HSV668 |
| 14.0 mm woven wire | HSV569 | HSV669 |
| 16.0mm woven wire | HSV570 | HSV670 |
| 18.0 mm woven wire | HSV571 | HSV671 |
| 19.0 mm woven wire | HSV572 | HSV672 |
| 20.0 mm woven wire | HSV573 | HSV673 |
| 22.4mm woven wire | HSV574 | HSV674 |
| 25.0 mm woven wire | HSV575 | HSV675 |
| 26.5 mm woven wire | HSV576 | HSV676 |
| 28.0mm woven wire | HSV577 | HSV677 |
| 31.5 mm woven wire | HSV578 | HSV678 |
| 35.5 mm woven wire | HSV579 | HSV679 |
| 37.5mm woven wire | HSV580 | HSV680 |
| 40.0 mm woven wire | HSV581 | HSV681 |
| 45.0mm woven wire | HSV582 | HSV682 |
| 50.0mm woven wire | HSV583 | HSV683 |
| 53.0mm woven wire | HSV584 | HSV684 |
| 56.0 mm woven wire | HSV585 | HSV685 |
| 63.0mm woven wire | HSV586 | HSV686 |
| 71.0 mm woven wire | HSV587 | HSV687 |
| 75.0mm woven wire | HSV588 | HSV688 |
| 80.0 mm woven wire | HSV589 | HSV689 |
| 90.0 mm woven wire | HSV590 | HSV690 |
| 100.0 mm woven wire | HSV591 | HSV691 |
| 106.0mm woven wire | HSV592 | HSV692 |
| 112.0 mm woven wire | HSV593 | HSV693 |
| 125.0mm woven wire | HSV594 | HSV694 |



European Standard Sieves to ISO 3310/2 and BS 410/2

Half Height Frames

Perforated Plate

| Aperture | – Frame Diameter – | | | |
|----------|--------------------|--------|------------|---------|
| | Square Hole | | Round Hole | |
| | 200mm | 300mm | 200mm | 300mm |
| 4.00 mm | HSV695 | HSV736 | HSV695R | HSV736R |
| 4.50 mm | HSV696 | HSV737 | HSV696R | HSV737R |
| 4.75 mm | HSV697 | HSV738 | HSV697R | HSV738R |
| 5.00 mm | HSV698 | HSV739 | HSV698R | HSV739R |
| 5.60 mm | HSV699 | HSV740 | HSV699R | HSV740R |
| 6.30 mm | HSV700 | HSV741 | HSV700R | HSV741R |
| 6.70 mm | HSV701 | HSV742 | HSV701R | HSV742R |
| 7.10mm | HSV702 | HSV743 | HSV702R | HSV743R |
| 8.00 mm | HSV703 | HSV744 | HSV703R | HSV744R |
| 9.00 mm | HSV704 | HSV745 | HSV704R | HSV745R |
| 9.50 mm | HSV705 | HSV746 | HSV705R | HSV746R |
| 10.0mm | HSV706 | HSV747 | HSV706R | HSV747R |
| 11.2 mm | HSV707 | HSV748 | HSV707R | HSV748R |
| 12.5 mm | HSV708 | HSV749 | HSV708R | HSV749R |
| 13.2 mm | HSV709 | HSV750 | HSV709R | HSV750R |
| 14.0 mm | HSV710 | HSV751 | HSV710R | HSV751R |
| 16.0mm | HSV711 | HSV752 | HSV711R | HSV752R |
| 18.0 mm | HSV712 | HSV753 | HSV712R | HSV753R |
| 19.0 mm | HSV713 | HSV754 | HSV713R | HSV754R |
| 20.0 mm | HSV714 | HSV755 | HSV714R | HSV755R |
| 22.4 mm | HSV715 | HSV756 | HSV715R | HSV756R |

| Aperture | – Frame Diameter – | | | |
|----------|--------------------|--------|------------|---------|
| | Square Hole | | Round Hole | |
| | 200mm | 300mm | 200mm | 300mm |
| 25.0 mm | HSV716 | HSV757 | HSV716R | HSV757R |
| 26.5 mm | HSV717 | HSV758 | HSV717R | HSV758R |
| 28.0 mm | HSV718 | HSV759 | HSV718R | HSV759R |
| 31.5 mm | HSV719 | HSV760 | HSV719R | HSV760R |
| 35.5 mm | HSV720 | HSV761 | HSV720R | HSV761R |
| 37.5 mm | HSV721 | HSV762 | HSV721R | HSV762R |
| 40.0 mm | HSV722 | HSV763 | HSV722R | HSV763R |
| 45.0 mm | HSV723 | HSV764 | HSV723R | HSV764R |
| 50.0 mm | HSV724 | HSV765 | HSV724R | HSV765R |
| 53.0 mm | HSV725 | HSV766 | HSV725R | HSV766R |
| 56.0 mm | HSV726 | HSV767 | HSV726R | HSV767R |
| 63.0 mm | HSV727 | HSV768 | HSV727R | HSV768R |
| 71.0 mm | HSV728 | HSV769 | HSV728R | HSV769R |
| 75.0 mm | HSV729 | HSV770 | HSV729R | HSV770R |
| 80.0 mm | HSV730 | HSV771 | HSV730R | HSV771R |
| 90.0 mm | HSV731 | HSV772 | HSV731R | HSV772R |
| 100mm | HSV732 | HSV773 | HSV732R | HSV773R |
| 106mm | HSV733 | HSV774 | HSV733R | HSV774R |
| 112mm | HSV734 | HSV775 | HSV734R | HSV775R |
| 125mm | HSV735 | HSV776 | HSV735R | HSV776R |



American Standard Sieves to ASTM E11-04

Woven Wire Mesh

| Aperture | – Frame Diameter – | |
|----------------------------|--------------------|---------|
| | 8 inch | 12 inch |
| 20µm - no. 635 woven wire | SV780 | SV857 |
| 25µm - no. 500 woven wire | SV781 | SV858 |
| 32µm - no. 450 woven wire | SV782 | SV859 |
| 38µm - no. 400 woven wire | SV783 | SV860 |
| 45µm - no. 325 woven wire | SV784 | SV861 |
| 53µm - no. 270 woven wire | SV785 | SV862 |
| 63µm - no. 230 woven wire | SV786 | SV863 |
| 75µm - no. 200 woven wire | SV787 | SV864 |
| 90µm - no. 170 woven wire | SV788 | SV865 |
| 106µm - no. 140 woven wire | SV789 | SV866 |
| 125µm - no. 120 woven wire | SV790 | SV867 |
| 150µm - no. 100 woven wire | SV791 | SV868 |
| 180µm - no. 80 woven wire | SV792 | SV869 |
| 212µm - no. 70 woven wire | SV793 | SV870 |
| 250µm - no. 60 woven wire | SV794 | SV871 |
| 300µm - no. 50 woven wire | SV795 | SV872 |
| 355µm - no. 45 woven wire | SV796 | SV873 |
| 425µm - no. 40 woven wire | SV797 | SV874 |
| 500µm - no. 35 woven wire | SV798 | SV875 |
| 600µm - no. 30 woven wire | SV799 | SV876 |
| 710µm - no. 25 woven wire | SV800 | SV877 |
| 850µm - no. 20 woven wire | SV801 | SV878 |
| 1.00mm - no. 18 woven wire | SV802 | SV879 |
| 1.18mm - no. 16 woven wire | SV803 | SV880 |
| 1.40mm - no. 14 woven wire | SV804 | SV881 |
| 1.70mm - no. 12 woven wire | SV805 | SV882 |
| 2.00mm - no. 10 woven wire | SV806 | SV883 |
| 2.36mm - no. 8 woven wire | SV807 | SV884 |

| Aperture | – Frame Diameter – | |
|--------------------------------|--------------------|---------|
| | 8 inch | 12 inch |
| 2.80mm - no. 7 woven wire | SV808 | SV885 |
| 3.35mm - no. 6 woven wire | SV809 | SV886 |
| 4.00mm - no. 5 woven wire | SV810 | SV887 |
| 4.75mm - no. 4 woven wire | SV811 | SV888 |
| 5.60mm - no. 3 1/2 woven wire | SV812 | SV889 |
| 6.30mm - 1/4 inch woven wire | SV813 | SV890 |
| 6.70mm - 0.265 inch woven wire | SV814 | SV891 |
| 8.00mm - 5/16 inch woven wire | SV815 | SV892 |
| 9.50mm - 3/8 inch woven wire | SV816 | SV893 |
| 11.2mm - 7/16 inch woven wire | SV817 | SV894 |
| 12.5mm - 1/2 inch woven wire | SV818 | SV895 |
| 13.2mm - 0.530 inch woven wire | SV819 | SV896 |
| 16.0mm - 5/8 inch woven wire | SV820 | SV897 |
| 19.0mm - 3/4 inch woven wire | SV821 | SV898 |
| 22.4mm - 7/8 inch woven wire | SV822 | SV899 |
| 25.0mm - 1 inch woven wire | SV823 | SV900 |
| 26.5mm - 1.06 inch woven wire | SV824 | SV901 |
| 31.5mm - 1 1/4 inch woven wire | SV825 | SV902 |
| 37.5mm - 1 1/2 inch woven wire | SV826 | SV903 |
| 45.0mm - 1 3/4 inch woven wire | SV827 | SV904 |
| 50.0mm - 2 inch woven wire | SV828 | SV905 |
| 53.0mm - 2.12 inch woven wire | SV829 | SV906 |
| 63.0mm - 2 1/2 inch woven wire | SV830 | SV907 |
| 75.0mm - 3 inch woven wire | SV831 | SV908 |
| 90.0mm - 3 1/2 inch woven wire | SV832 | SV909 |
| 100mm - 4 inch woven wire | SV833 | SV910 |
| 106mm - 4.24 inch woven wire | SV834 | SV911 |
| 125mm - 5 inch woven wire | SV835 | SV912 |

American Standard Sieves to ASTM E11-04

Half Height Frames

Woven Wire Mesh

| Aperture | – Frame Diameter – | |
|----------------------------|--------------------|---------|
| | 8 inch | 12 inch |
| 20µm - no. 635 woven wire | HSV780 | HSV857 |
| 25µm - no. 500 woven wire | HSV781 | HSV858 |
| 32µm - no. 450 woven wire | HSV782 | HSV859 |
| 38µm - no. 400 woven wire | HSV783 | HSV860 |
| 45µm - no. 325 woven wire | HSV784 | HSV861 |
| 53µm - no. 270 woven wire | HSV785 | HSV862 |
| 63µm - no. 230 woven wire | HSV786 | HSV863 |
| 75µm - no. 200 woven wire | HSV787 | HSV864 |
| 90µm - no. 170 woven wire | HSV788 | HSV865 |
| 106µm - no. 140 woven wire | HSV789 | HSV866 |
| 125µm - no. 120 woven wire | HSV790 | HSV867 |
| 150µm - no. 100 woven wire | HSV791 | HSV868 |
| 180µm - no. 80 woven wire | HSV792 | HSV869 |
| 212µm - no. 70 woven wire | HSV793 | HSV870 |
| 250µm - no. 60 woven wire | HSV794 | HSV871 |
| 300µm - no. 50 woven wire | HSV795 | HSV872 |
| 355µm - no. 45 woven wire | HSV796 | HSV873 |
| 425µm - no. 40 woven wire | HSV797 | HSV874 |
| 500µm - no. 35 woven wire | HSV798 | HSV875 |
| 600µm - no. 30 woven wire | HSV799 | HSV876 |
| 710µm - no. 25 woven wire | HSV800 | HSV877 |
| 850µm - no. 20 woven wire | HSV801 | HSV878 |
| 1.00mm - no. 18 woven wire | HSV802 | HSV879 |
| 1.18mm - no. 16 woven wire | HSV803 | HSV880 |
| 1.40mm - no. 14 woven wire | HSV804 | HSV881 |
| 1.70mm - no. 12 woven wire | HSV805 | HSV882 |
| 2.00mm - no. 10 woven wire | HSV806 | HSV883 |
| 2.36mm - no. 8 woven wire | HSV807 | HSV884 |

| Aperture | – Frame Diameter – | |
|--------------------------------|--------------------|---------|
| | 8 inch | 12 inch |
| 2.80mm - no. 7 woven wire | HSV808 | HSV885 |
| 3.35mm - no. 6 woven wire | HSV809 | HSV886 |
| 4.00mm - no. 5 woven wire | HSV810 | HSV887 |
| 4.75mm - no. 4 woven wire | HSV811 | HSV888 |
| 5.60mm - no. 3 1/2 woven wire | HSV812 | HSV889 |
| 6.30mm - 1/4 inch woven wire | HSV813 | HSV890 |
| 6.70mm - 0.265 inch woven wire | HSV814 | HSV891 |
| 8.00mm - 5/16 inch woven wire | HSV815 | HSV892 |
| 9.50mm - 3/8 inch woven wire | HSV816 | HSV893 |
| 11.2mm - 7/16 inch woven wire | HSV817 | HSV894 |
| 12.5mm - 1/2 inch woven wire | HSV818 | HSV895 |
| 13.2mm - 0.530 inch woven wire | HSV819 | HSV896 |
| 16.0mm - 5/8 inch woven wire | HSV820 | HSV897 |
| 19.0mm - 3/4 inch woven wire | HSV821 | HSV898 |
| 22.4mm - 7/8 inch woven wire | HSV822 | HSV899 |
| 25.0mm - 1 inch woven wire | HSV823 | HSV900 |
| 26.5mm - 1.06 inch woven wire | HSV824 | HSV901 |
| 31.5mm - 1 1/4 inch woven wire | HSV825 | HSV902 |
| 37.5mm - 1 1/2 inch woven wire | HSV826 | HSV903 |
| 45.0mm - 1 3/4 inch woven wire | HSV827 | HSV904 |
| 50.0mm - 2 inch woven wire | HSV828 | HSV905 |
| 53.0mm - 2.12 inch woven wire | HSV829 | HSV906 |
| 63.0mm - 2 1/2 inch woven wire | HSV830 | HSV907 |
| 75.0mm - 3 inch woven wire | HSV831 | HSV908 |
| 90.0mm - 3 1/2 inch woven wire | HSV832 | HSV909 |
| 100mm - 4 inch woven wire | HSV833 | HSV910 |
| 106mm - 4.24 inch woven wire | HSV834 | HSV911 |
| 125mm - 5 inch woven wire | HSV835 | HSV912 |

Lids and Receivers

| Lid | – Frame Diameter – | |
|----------------------------|--------------------|---------|
| | 8 inch | 12 inch |
| Lid | SV930 | SV935 |
| Receiver | SV931 | SV936 |
| Receiver/Separator | SV933 | SV937 |
| Wet Washing Lid & Receiver | SV458 | SV459 |



Grain Sieves to ISO 5223

Full Height Frames

| Aperture | – Frame Diameter – | |
|-------------------------------|--------------------|--|
| | 200mm | |
| 1.0mm x 20.0mm slotted plate | SG001 | |
| 1.7mm x 20.0mm slotted plate | SG002 | |
| 1.8mm x 20.0mm slotted plate | SG003 | |
| 1.9mm x 20.0mm slotted plate | SG004 | |
| 2.0mm x 20.0mm slotted plate | SG005 | |
| 2.20mm x 20.0mm slotted plate | SG006 | |
| 2.50mm x 20.0mm slotted plate | SG007 | |
| 2.80mm x 20.0mm slotted plate | SG008 | |
| 3.55mm x 20.0mm slotted plate | SG009 | |
| 4.50mm round hole | SG010 | |

Half Height Frames

| Aperture | – Frame Diameter – | |
|-------------------------------|--------------------|--|
| | 200mm | |
| 1.00mm x 20.0mm slotted plate | HSG001 | |
| 1.70mm x 20.0mm slotted plate | HSG002 | |
| 1.80mm x 20.0mm slotted plate | HSG003 | |
| 1.90mm x 20.0mm slotted plate | HSG004 | |
| 2.00mm x 20.0mm slotted plate | HSG005 | |
| 2.20mm x 20.0mm slotted plate | HSG006 | |
| 2.50mm x 20.0mm slotted plate | HSG007 | |
| 2.80mm x 20.0mm slotted plate | HSG008 | |
| 3.55mm x 20.0mm slotted plate | HSG009 | |
| 4.50mm round hole | HSG010 | |

Large Capacity Screen Tray Sieves - 457 x 660 x 75mm

Woven Wire Mesh

Aperture

| | |
|------------------|-------|
| 20µm woven wire | SZ001 |
| 25µm woven wire | SZ002 |
| 32µm woven wire | SZ003 |
| 36µm woven wire | SZ004 |
| 38µm woven wire | SZ005 |
| 40µm woven wire | SZ006 |
| 45µm woven wire | SZ007 |
| 50µm woven wire | SZ008 |
| 53µm woven wire | SZ009 |
| 56µm woven wire | SZ010 |
| 63µm woven wire | SZ011 |
| 71µm woven wire | SZ012 |
| 75µm woven wire | SZ013 |
| 80µm woven wire | SZ014 |
| 90µm woven wire | SZ015 |
| 100µm woven wire | SZ016 |
| 106µm woven wire | SZ017 |
| 112µm woven wire | SZ018 |
| 125µm woven wire | SZ019 |
| 140µm woven wire | SZ020 |
| 150µm woven wire | SZ021 |
| 160µm woven wire | SZ022 |
| 180µm woven wire | SZ023 |
| 200µm woven wire | SZ024 |
| 212µm woven wire | SZ025 |
| 224µm woven wire | SZ026 |
| 250µm woven wire | SZ027 |
| 280µm woven wire | SZ028 |
| 300µm woven wire | SZ029 |
| 315µm woven wire | SZ030 |
| 355µm woven wire | SZ031 |
| 400µm woven wire | SZ032 |
| 425µm woven wire | SZ033 |

Aperture

| | |
|-------------------|-------|
| 450µm woven wire | SZ034 |
| 500µm woven wire | SZ035 |
| 560µm woven wire | SZ036 |
| 600µm woven wire | SZ037 |
| 630µm woven wire | SZ038 |
| 710µm woven wire | SZ039 |
| 800µm woven wire | SZ040 |
| 850µm woven wire | SZ041 |
| 900µm woven wire | SZ042 |
| 1.00mm woven wire | SZ043 |
| 1.12mm woven wire | SZ044 |
| 1.18mm woven wire | SZ045 |
| 1.25mm woven wire | SZ046 |
| 1.40mm woven wire | SZ047 |
| 1.60mm woven wire | SZ048 |
| 1.70mm woven wire | SZ049 |
| 1.80mm woven wire | SZ050 |
| 2.00mm woven wire | SZ051 |
| 2.24mm woven wire | SZ052 |
| 2.36mm woven wire | SZ053 |
| 2.50mm woven wire | SZ054 |
| 2.80mm woven wire | SZ055 |
| 3.15mm woven wire | SZ056 |
| 3.35mm woven wire | SZ057 |
| 3.55mm woven wire | SZ058 |
| 4.00mm woven wire | SZ059 |
| 4.50mm woven wire | SZ060 |
| 4.75mm woven wire | SZ061 |
| 5.00mm woven wire | SZ062 |
| 5.60mm woven wire | SZ063 |
| 6.30mm woven wire | SZ064 |
| 6.70mm woven wire | SZ065 |
| 7.10mm woven wire | SZ066 |

Aperture

| | |
|-------------------|-------|
| 8.00mm woven wire | SZ067 |
| 9.00mm woven wire | SZ068 |
| 9.50mm woven wire | SZ069 |
| 10.0mm woven wire | SZ070 |
| 11.2mm woven wire | SZ071 |
| 12.5mm woven wire | SZ072 |
| 13.2mm woven wire | SZ073 |
| 14.0mm woven wire | SZ074 |
| 16.0mm woven wire | SZ075 |
| 18.0mm woven wire | SZ076 |
| 19.0mm woven wire | SZ077 |
| 20.0mm woven wire | SZ078 |
| 22.4mm woven wire | SZ079 |
| 25.0mm woven wire | SZ080 |
| 26.5mm woven wire | SZ081 |
| 28.0mm woven wire | SZ082 |
| 31.5mm woven wire | SZ083 |
| 35.5mm woven wire | SZ084 |
| 37.5mm woven wire | SZ085 |
| 40.0mm woven wire | SZ086 |
| 45.0mm woven wire | SZ087 |
| 50.0mm woven wire | SZ088 |
| 53.0mm woven wire | SZ089 |
| 56.0mm woven wire | SZ090 |
| 63.0mm woven wire | SZ091 |
| 71.0mm woven wire | SZ092 |
| 75.0mm woven wire | SZ093 |
| 80.0mm woven wire | SZ094 |
| 90.0mm woven wire | SZ095 |
| 100mm woven wire | SZ096 |
| 106mm woven wire | SZ097 |
| 112mm woven wire | SZ098 |
| 125mm woven wire | SZ099 |

Dustpan

SZ141

Perforated Plate

Aperture

| | |
|-------------------------|-------|
| 4.00mm perforated plate | SZ100 |
| 4.50mm perforated plate | SZ101 |
| 4.75mm perforated plate | SZ102 |
| 5.00mm perforated plate | SZ103 |
| 5.60mm perforated plate | SZ104 |
| 6.30mm perforated plate | SZ105 |
| 6.70mm perforated plate | SZ106 |
| 7.10mm perforated plate | SZ107 |
| 8.00mm perforated plate | SZ108 |
| 9.00mm perforated plate | SZ109 |
| 9.50mm perforated plate | SZ110 |
| 10.0mm perforated plate | SZ111 |
| 11.2mm perforated plate | SZ112 |
| 12.5mm perforated plate | SZ113 |

Aperture

| | |
|-------------------------|-------|
| 13.2mm perforated plate | SZ114 |
| 14.0mm perforated plate | SZ115 |
| 16.0mm perforated plate | SZ116 |
| 18.0mm perforated plate | SZ117 |
| 19.0mm perforated plate | SZ118 |
| 20.0mm perforated plate | SZ119 |
| 22.4mm perforated plate | SZ120 |
| 25.0mm perforated plate | SZ121 |
| 26.5mm perforated plate | SZ122 |
| 28.0mm perforated plate | SZ123 |
| 31.5mm perforated plate | SZ124 |
| 35.5mm perforated plate | SZ125 |
| 37.5mm perforated plate | SZ126 |
| 40.0mm perforated plate | SZ127 |

Aperture

| | |
|-------------------------|-------|
| 45.0mm perforated plate | SZ128 |
| 50.0mm perforated plate | SZ129 |
| 53.0mm perforated plate | SZ130 |
| 56.0mm perforated plate | SZ131 |
| 63.0mm perforated plate | SZ132 |
| 71.0mm perforated plate | SZ133 |
| 75.0mm perforated plate | SZ134 |
| 80.0mm perforated plate | SZ135 |
| 90.0mm perforated plate | SZ136 |
| 100mm perforated plate | SZ137 |
| 106mm perforated plate | SZ138 |
| 112mm perforated plate | SZ139 |
| 125mm perforated plate | SZ140 |



Wet Washing Sieves to ISO 3310/1 and ASTM E11-04

200mm diameter

| Aperture | Depth | |
|-----------------|-------|-------|
| 63µm woven wire | 200mm | SV444 |
| 75µm woven wire | 200mm | SV445 |
| 63µm woven wire | 100mm | SV443 |
| 75µm woven wire | 100mm | SV446 |

450mm diameter

| Aperture | Depth | |
|------------------|-------|-------|
| 63µm woven wire | 300mm | SV460 |
| 75µm woven wire | 300mm | SV455 |
| 150µm woven wire | 300mm | SV450 |

300mm diameter

| Aperture | Depth | |
|-----------------|-------|-------|
| 63µm woven wire | 300mm | SV448 |
| 75µm woven wire | 300mm | SV447 |
| 63µm woven wire | 225mm | SV449 |
| 75µm woven wire | 225mm | SV494 |

8 inch diameter

| Aperture | Depth | |
|-----------------|----------|-------|
| 63µm woven wire | 4 inches | SV917 |
| 75µm woven wire | 4 inches | SV918 |
| 63µm woven wire | 8 inches | SV919 |
| 75µm woven wire | 8 inches | SV920 |

Test Sieve Specifications

Frame Diameter

Frame Material

Sieve Depth

| | | |
|--------------------|--|-----------|
| 75mm | Stainless Steel | 40mm |
| 100mm | Stainless Steel | 40mm |
| 150mm | Stainless Steel | 50mm |
| 200mm | Stainless Steel | 50mm |
| 250mm | Stainless Steel | 50mm |
| 300mm | Stainless Steel | 75mm |
| 315mm | Stainless Steel | 75mm |
| 350mm | Stainless Steel | 65mm |
| 400mm | Stainless Steel | 75mm |
| 450mm | Plated Steel | 100mm |
| 8 inch | Stainless Steel | 2 inches |
| 12 inch | Stainless Steel | 3 inches |
| Wet Washing Sieves | 200mm, 300mm, 8 and 12 inch diameter - Stainless Steel | see above |
| Wet Washing Sieves | 450mm diameter Plated Steel | 300mm |
| 200mm Grain Sieve | Stainless Steel | 50mm |

The sieve diameter is the length of the straight line from inside edge to inside edge of the frame, passing through the centre of the sieve.

All woven wire mesh used in Impact Test Sieves is manufactured from stainless steel.

The depth of the sieve is measured from the top of the frame to the surface of the woven wire mesh or perforated plate.

Perforated plates used in Impact Test Sieves are manufactured from mild steel. The only exception to this is our range of Grain Sieves, manufactured with stainless steel perforated plates