

## BALTEC Professional Analytical Balance (0.1mg)

### ■ Application Mode

- ✧ Percentage Weighing
- ✧ With density determination
- ✧ With dynamic weighing function
- ✧ Selectable weighing and count mode
- ✧ Pricing function
- ✧ High-low weighing function
- ✧ Gross/Net/Tare weight weighing
- ✧ Peak holding function
- ✧ Accumulate function


Weighing functions are selectable on PC. Measurements can be input into EXCEL, WORD and other files on PC.

### ■ Characteristic

- ◆ Super Monobloc sensor
- ◆ Internal calibration with one key touch.
- ◆ Real-time temperature compensation
- ◆ Seven level anti-vibration adjustable
- ◆ VATN Backlit LCD Display
- ◆ Weighing stableness time adjustable
- ◆ More than 20unit shield: g, mg, ct, kg...
- ◆ Default weighing mode automatically
- ◆ Date, time and temperature setting
- ◆ With carrying handle and rear tool box
- ◆ With underhook for below weighing
- ◆ Measurements are directly readable on PC
- ◆ Anti-static glass draftshield
- ◆ Overload protection
- ◆ RS232 and USB port can communicate with external equipment

Our sensor technology's features are fast weighing response, high durability, continually coherence, compatible of sensitivity and stability.

### Internal CAL With One-key Touch



Item No.	Capacity	Readability	Repeatability	Linearity	Pan Size(mm)	Response Time	Housing(mm)
BM204-I	220g	0.1mg	±0.1mg	±0.2mg	Φ90	2~3 secs	345×223×331

\* All Specifications & Dimensions are subject to change without Prior notice.