

# TB60 Single Punch Tablet Machine

## MILLMIX



TB60 worm gear single punch tablet machine is a small single punch tablet machine. The machine has good performance, easy to operate and easy to move. Various granular, crystalline or powdery raw materials with good fluidity can be compressed into tablets, which are suitable for small batch production and teaching use in scientific research experiments. Time-saving and labor-saving during electric operation, high efficiency; this machine is only equipped with one set of punching mold, the filling depth and pressing thickness of the material can be adjusted. (adjusted according to the weight, thickness and hardness of the product) different products can be pressed simply by changing different molds. It is a relatively lightweight electric tablet press.

### Main Technical Parameters

Model	TB60
Tableting Pressure	0-60 KN
Diameter of Pressing Tablet	20 mm
Max. Filling Depth	18 mm
Max. Thickness of Pressing Tablet	7 mm
Capacity	4000 pcs/h
Power	1.5 kw
Dimension	500x410x680 mm
Weight	85 kg
Voltage	220VAC / 50Hz