

BALTEC BALANCE Large Capacity Precision Balance 0.01g

Application Mode

- ✘ Percentage Weighing
- ✘ With density determination
- ✘ With dynamic weighing function
- ✘ Selectable weighing and count mode
- ✘ Pricing function
- ✘ High-low weighing function
- ✘ Gross/Net/Tare weight weighing
- ✘ Peak holding function
- ✘ Accumulate function

Characteristic

- ◆ Loadcell Sensor
- ◆ Mostly Auto dual range/dual precision
- ◆ External Calibration
- ◆ Real-time temperature compensation
- ◆ Seven level anti-vibration adjustable
- ◆ Anti-static mental shell
- ◆ VATN Backlit LCD Display
- ◆ Weighing stableness time adjustable
- ◆ More than 20unit shield: g, kg, ct...
- ◆ Default weighing mode automatically
- ◆ Date, time and temperature setting
- ◆ With carrying handle
- ◆ RS232 to communicate with PC and printer
- ◆ Measurements are directly readable on PC
- ◆ Optional with underhook
- ◆ Optional with rechargeable battery

External CAL (With CAL weights) Loadcell Sensor

BL serials are excellent products which are made by higher level technology of 0.001g-0.1g accuracy. It has firmly iron shell, transparent and beautiful glass windproof which is not used in same level product. It has quick response, high sensitivity and stability.



Characteristic



Application Mode



Item No.	Capacity	Readability	Repeatability	Linearity	Pan Size(mm)	Response Time	Housing(mm)
BL3002	3200g/4200 g	0.01 g/0.0 5g	±0.02 g/±0.05 g	±0.03 g/±0.1 g	168 x 168	2-3 Sec.	292 x 208 x 88
BL4002	4000g	0.01g	±0.02g	±0.03g	168 x 168	2-3 Sec.	292 x 208 x 88
BL5002	5000g	0.01g	±0.02g	±0.03g	168 x 168	2-3 Sec.	292 x 208 x 88
BL6002	6000g	0.01g	±0.02g	±0.03g	168 x 168	2-3 Sec.	292 x 208 x 88
BL8002	8000g	0.01g	±0.02g	±0.03g	168 x 168	2-3 Sec.	292 x 208 x 88