

Product name	Chiller
Model	CB10-40
Volume (L)	10
Temperature range (°C)	-40 to Room temperature
Temperature stability (°C)	0.5
Temperature resolution (°C)	0.1
Display	LED
Safety protection	Delay, phase sequence, over-current, overheat
Refrigerant	R404A
Compressor power	1.5HP
Refrigerate capacity (W)	2548-980W
Pump flow rate (L/min)	32L/min
Open bath diameter (mm)	290
Chamber	SUS304
Overall size (mm)	710 x 660 x 890
Power voltage	Single phase 220V, 50Hz

