

BALTEC BALANCE

BC Precision Scale (0.01g)

■ Application Mode

- ✧ Percentage Weighing
- ✧ Selectable weighing and count mode

Optional with RS232: Weighing functions are selectable on PC. Measurements can be input into EXCEL, WORD and other files on PC.

■ Characteristic

- ◆ Loadcell Sensor
- ◆ External Calibration
- ◆ Six level anti-vibration adjustable
- ◆ Weighing stableness time adjustable
- ◆ More than 14unit shield: g, kg, ct...
- ◆ Default weighing mode automatically
- ◆ Optional with rechargeable battery
- ◆ Optional with RS232 to communicate with external devices, like PC and printer
- ◆ Optional with rechargeable battery
- ◆ Optional with external display
- ◆ Optional with underhook

External CAL (With CAL weights) Loadcell Sensor

BC is a standard product with the function of AC-DC which is convenient for outdoor users. It is LCD display and is optional with external display. It is the most cost-effective balances.



CAL




SENSOR



Characteristic



Application Mode



Model	Capacity(g)	Readability	Repeatability	Linearity(g)	Pan(mm)	Response Time	Housing(mm)
BC1002-E	1000g	0.01g	±0.01g	±0.02g	Φ130	3 sec.	295×208×88