

Characteristics :

- Fully stainless steel structure.
- Hand wheel type of quick-open door structure.
- Door safety lock system.
- Digital display of working status, touch of key.
- Auto discharge the cool air, and steam discharging automatically after sterilization.
- Over temperature & over pressure auto-protection.
- Safe protection of water lacking.
- Self-inflating type seal.



Model	CA-35	CA-50	CA-75
Volume	35L	50L	75L
Dimension	Ø318 x H450 mm	Ø340 x H550 mm	Ø400 x H600 mm
Working temperature	105-134°C		
Max working pressure	0.23 MPa		
Accuracy	±1°C		
Timer	0-99 Hours 59 Min		
Heater	2.5Kw	3Kw	4.5Kw
Power	AC220V / 50Hz	AC220V / 50Hz	AC220V / 50Hz
Overall Dimension	450x450x1010 mm.	510x470x1130 mm.	560x560x1120 mm.
Net Weight	56 Kg	68 Kg	80 Kg