www.shinjineng.com



The user can program up to 5 programmable preset cleaning cycles consisting of a time and temperature. As different tools and applications have different sterilization needs, our product lets you configure up to 5 presets for easy use.



Fully Automated Vertical Steam Sterilizer

5 Programmable Preset Slots

Chamber



Accessories



Vertical stainless steel wire baskets containers in different sizes for all autoclave models





Specifications

	<i>SJ-AD60</i>	<i>SJ-AD80</i>	<i>SJ-AD100</i>
Chamber Capacity	62Liter	81Liter	106Liter
Chamber Size (Ø x H)mm	350Ф x 650	400Ф x 650	400Φ x 850
Dimension(WxDxH)mm	700 x 750 x 825	700 x 800 x 825	700 x 800 x 1025
Heater Capacity	3Kw	4Kw	4Kw
Temp Range	121°C ~ 135°C		
Controller	LCD Microprocessor Programble control		
Chamber Material	Stainless steel (SUS 304)		
Max Pressure	0.263MPA		
Electronic Power	AC 220V, 50-60Hz Single-phase		
	Water level sensor / Safety Door interlock/		
Safety Device	Over-pressure prevention/ Over-heat prevention/		
	Open temp. sensor detection/ Safety valve		

Convenience



*All sizes can be custom made

* Nipple is made of stainless steel and does not rust.

LEVER

SJ-AD (lever type) sterilizer is convenient to use and offers easy loading and unloading with its moving lever mechanism.