

Table type low speed centrifuge

HIQ LAB



FEATURES :

- ◆ With brushless motor, maintenance-free.
- ◆ Controlled by Microprocessor with touching panel and digital display, easy to operate.
- ◆ Steel body, stainless steel centrifugal chamber, resist corrosion, easy to clean, safe and stable.
- ◆ Arbitrarily setting the parameter of time and speed etc, easy to use.
- ◆ 0-9 Accel/Decel modes, 40 self-defining programs can be stored and recall.
- ◆ Short spin function.
- ◆ It's mainly used in biochemistry Laboratory and Clinical Laboratory in Hospital and so on.

TECHNICAL PARAMETER :

Model	CENTRI-6300M
Max Speed	6000r/min
Max RCF	5120×g
Capacity	300ml
Speed Accuracy	±20r/min
Time Rang	0-999min/sec
Power Supply	AC220V 50Hz / AC110V 60Hz,15A
Dimension	340×400×270mm
Weight	20kg

ROTOR OPTIONS :

Rotor Type	Rotor No.	Capacity	Max Speed (r/min)	Max RCF (×g)	Dimension of tube (mm)
Swing Rotor Microplate Rotor	No.38	4×50ml	6000	5120	Φ28×114
	No.38	6×50ml	6000	5120	Φ28×114
	No.37	12×10ml	6000	5120	Φ16×82
	No.37	12×15ml	6000	5120	Φ16×115