

Orbital Shaker—Double Platform

Features

- Touch-screen panel. LED 2-screen digital display
- Actual and setting speed
- Shaking smoothly and quietly
- Suitable for carrying heavy duty
- With variable speed, timer control and continuous
- Operation
- Memory of power cut ,uses brushless AC motor
- Differential clamps available
- Double platform for large capacity



Specifications

Model No.	SH-202D	SH-602D	SH-902D
Control Panel	Digital		
Shaking Speed	20~150 r.p.m		
Shaking Width	25 mm		
Timer	999 h / 999 min / 999 sec		
Shaking Plate	400 x 300 2 pcs	580 x 460 2 pcs	960 x 600 2 pcs
Outside Dimension (WxDxH)	400 x 300 x 500 mm	600 x 480 x 500 mm	960 x 600 x 500 mm
Weight	20 kg	60 kg	107 kg

Capacity for flasks

Flask Clamp	125ml	250ml	500ml
SH-202D	40	24	12
SH-602D	96	60	40
SH-902D	192	120	80

Power Voltage : 110V/220V 50Hz/60Hz